



PRESSURE FIT SYSTEM—2 POST

Item #	Part No.	Description	Qty	✓ = Included O = Optional
1	433101	Pressure Fit Track & Trolley Assembly	1 piece	✓
2	433102	Pressure Fit Post Assembly	2 pieces	✓
3	433103	Pressure Fit Bottom Plate Assembly	2 pieces	✓
4	433104	Pressure Fit Top Plate Assembly	2 pieces	✓
5	633545	Pressure Fit Track End Cover	2 pieces	✓
6	341501	Pressure Fit (2 Post) Carry Bag	1 piece	O



- Built in Post Leveling system for ease of alignment.



- Pressure gauge indicator to help ensure proper installation.



- Built in Trolley slides smoothly

- Built in Trolley slides smoothly along the extended track.

