

Free-standing overhead systems

Waverley Glen **free-standing overhead systems** are everything free-standing systems were meant to be:

- they stand freely on the floor without support from floor-to-ceiling pressure posts
- they're collision resistant and don't require adjusting over time
- they're independent of ceiling type and wall type

The systems suit rental requirements as well as situations where ceiling or wall mounting is not possible or desirable. All free-standing systems are ideal for use with the Waverley Glen *transportable* or Griffin lift. Units are assembled with standard wrenches and Allan keys.



T-shape 3-post bedroom system

Components

- 1 transfer rail with trolley
- 1 repositioning rail with trolley
- 2 cross braces
- 2 bed posts with leveling feet
- 1 room post with beveled foot

Specifications

The system is 7 feet high without the leveling feet.

The repositioning track fits a bed length of up to 7 ft 3 in. The transfer track measures 4 ft 8 in.

Each post weighs 23.5 lbs. The tracks weigh 12 lbs and 8 lbs.

Waverley Glen Systems Ltd.

Canada | 87 Sharer Road, Vaughan ON L4L 8Z3 | T 1.800.265.0677 | 905.850.0093 | F 905.850.8377

USA | 2351 Grissom Drive, St. Louis, MO 63146 | T 1.866.891.6502 | 314.692.9135 | F 314.692.7858 | www.waverleyglen.com

Square 4-post area covering system

Components

- 2 support rails
- 2 support rails with trolleys
- 1 gantry rail
- 4 posts with leveling feet

Specifications

The system is 7 feet high without the leveling feet.

The gantry is a square measuring 8 ft 9 in along the outside and 8 ft 2 in along the inside.

Each track weighs 14 lbs and each post 23.5 lbs.



Cross-shape 3-post bathroom system

Components

- 1 transfer rail with trolley
- 1 support rail
- 2 bath posts with leveling feet
- 1 weighted room post with leveling feet counterweights

Specifications

The system is 7 feet high without the leveling feet.

The support track fits a minimum bath length of 4 ft 10 in. The transfer track measures 7 ft 3 in. The minimum distance between the weighted post and the support track is 5 ft.

Each track weighs 10.5 lbs, each bath post 23.5 lbs, the weighted post 53.5 lbs and each counterweight 7.3 lbs.

