


DISPOSABLE SLING

Waverley Glen is pleased to be able to offer disposable slings. With Infection Control becoming an area of greater focus and attention, Waverley Glen has the solution. Disposable Slings are designed to be used by one individual and disposed of when; a) the patient/resident no longer requires or b) the sling is soiled.

The disposable slings are based on the design, safety, comfort and application of the Waverley Glen Universal sling pattern.

All disposable slings are rated to 600 lbs., providing the same weight rating as all "standard" Waverley Glen slings.

Label indicates when the sling is in need of disposal. Typical safe use ranges from 1-3 months (dependent upon frequency of use)



CUSTOM SLING

Although we offer over 100 standard slings, everyone is unique and their needs are different. That is why we will customize any of our standard slings to help suit your needs better.

You can modify an existing sling design, like adding an extra strap, or we can start right from scratch.

Before creating a custom sling, our engineering department reviews the design to ensure that there will be no loss of integrity with the changes.

Once the custom sling has been made by our in-house sling manufacturing, we weight test each custom sling to ensure its strength.

Common customizations include adding: loops, straps and padding.

WE'RE THE COMPANY THAT BRINGS YOU LYNETTE...



AND LYNETTE BRINGS YOU THE WAVERLEY GLEN SLING SELECTION

Waverley Glen slings are recognized by health care professionals throughout North America for exceptional fit and durability. Our range of sling styles and sizes will accommodate virtually any lift, positioning and transfer requirement. And, if our standard selection does not provide a sling solution, we will custom design a sling system to your specifications.

Waverley Glen slings are manufactured from durable nylon mesh, polyester fabrics and quilted padding to provide contoured support, comfort and safety. As well, these long-wearing, quick-drying slings contain no metal parts, making them convenient to launder and easy to return to service, without the need to locate mislaid inserts.

For more information visit www.waverleyglen.com

WAVERLEY GLEN
A PRISM MEDICAL COMPANY

Waverley Glen Systems Ltd.

87 Sharer Road, Vaughan, ON Canada L4L 8Z3 | T 905.660.0213 | 1.800.265.0677 | www.waverleyglen.com


WAVERLEY GLEN
A PRISM MEDICAL COMPANY

WE'RE THE COMPANY THAT BRINGS YOU A WIDE SELECTION OF SLINGS




HAMMOCK SLING

Acknowledged by professionals as the most comfortable and supportive sling for both home and institutional use, the versatile Hammock Sling combines toileting access with full body and head support. The sling comes with leg and head straps and features loops to allow for various sitting positions. The Hammock is available in either regular mesh or with quilted material.




HYGIENE SLING

The Hygiene Sling was specifically designed for toileting and hygiene functions. Access is unrestricted for most users, making this sling ideal for independent transfer and toileting. The Hygiene Sling features quilted material and comes with leg and shoulder straps with loops for various seating positions. A wide velcro strap across the chest provides additional upper body support.




UNIVERSAL SLING

The Universal sling is easy-to-use and easy-to-fit. It is a multi-purpose sling that provides toileting access as well as good back and thigh support. The sling comes with leg, hip and shoulder straps and features loops to accommodate various seating positions. Available in either regular mesh or with quilted material, the Universal Sling can also be supplied with added head and neck support.




STAND-AID SLING

Designed by Waverley Glen's expert research team, the Stand-Aid sling will comfortably wrap around the patients waist and with limited pressure will help raise the person to a standing position. The Stand-Aid sling is padded for extra comfort and comes with Waverley Glen's signature 3 loop design, which allows for greater options for comfort when performing the lift.



POSITIONING SLING

The Positioning Sling is intended to be permanently in place as part of the bedding. The sling has four pairs of straps and features loops to support the head in various positions. Re-positioning and turning, done manually under time pressure, poses a significant risk of injury to caregivers. The Positioning Sling is designed to address that risk. The sling is also used for transfers to other horizontal surfaces.



WALKING/STANDING SLING

The Walking sling was developed with the help of physiotherapists, this supportive sling assists in the safe lifting of individuals from a seated to a standing position. The Walking Sling offers upper body support and freedom of movement to individuals performing exercises in a rehabilitation setting.

The Standing sling is a dual purpose sling affording a natural sit to stand movement. Functional from both bedside and chair surface, this sling encourages pivot transfers, ambulation and toileting.

